

Product Datasheet

Mouse LIF Protein (orb1146996)

Description	Mouse LIF Protein, Tag Free is expressed from human 293 cells (HEK293). It contains AA Ser 24 - Phe ...
Reactivity	Mouse
Endotoxins	0.1 EU per µg
Conjugation	Unconjugated
Target	LIF
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 20.0 kDa. The protein migrates as 35-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Mouse LIF Protein, Tag Free (orb1146996) is expressed from human 293 cells (HEK293). It contains AA Ser 24 - Phe 203 (Accession # P09056-1).
MW	20.0 kDa
Uniprot ID	P09056
NCBI	NP_032527.1
Expiration Date	6 months from date of receipt.