

Product Datasheet

Mouse CD36 / SR-B3 Protein (orb1147017)

Description	Mouse CD36, His Tag is expressed from human 293 cells (HEK293). It contains AA Gly 30 - Lys 439 (Acc...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD36
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 48.3 kDa. The protein migrates as 66-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_031669.3
Purity	95%
Source	Mouse CD36, His Tag (orb1147017) is expressed from human 293 cells (HEK293). It contains AA Gly 30 - Lys 439 (Accession # Q08857-1).
MW	48.3 kDa
Uniprot ID	Q08857
NCBI	NP_031669.3
Expiration Date	6 months from date of receipt.