

Product Datasheet

Z-VAD-FMK (orb1147018)

Description	Pan-caspase inhibitor; Peptides.
Form/Appearance	White lyophilised solid
Storage	Store dry, dark and frozen
Note	For research use only
Protein Sequence	Z-VAD-FMK
Formula	$C_{22}H_{30}FN_3O_7$
Purity	> 98% by TLC
MW	467.49 Da
CAS Number	187389-52-2
Solubility (25°C)	Soluble to 10 mg/mL in DMSO
Expiration Date	12 months from date of receipt.