

Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleocapsid Protein (orb1147171)

Description	Human SARS-COV2 Nucleocapsid protein Recombinantly expressed in mammalian expression system.
Reactivity	Human
Preservatives	PBS
Form/Appearance	Lyophilized/ Liquid
Storage	Store in -20°C
Tag	C-terminus 6× Histag
Note	For research use only
Application notes	Covid related research
Protein Sequence	Full length
Purity	85%, SDS-PAGE (Coomassie blue staining)
Source	E.coli
MW	51kDa
Expiration Date	6 months from date of receipt.