

Product Datasheet

SARS-CoV-2 (2019-nCoV) S Protien full (orb1147172)

Description	Human SARS-COV2 S protein Recombinantly expressed in mammalian expression system.
Reactivity	Human
Preservatives	PBS
Form/Appearance	Lyophilized/ Liquid
Storage	Store in -20°C
Tag	C-terminus 6× Histag
Note	For research use only
Application notes	Covid related research
Protein Sequence	Full length
Purity	90%, SDS-PAGE (Coomassie blue staining)
Source	E.coli
MW	138KDa
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.