

Product Datasheet

SARS-CoV-2 (2019-nCoV) Spike S1 protein (orb1147188)

Description	Human SARS-COV2 S1 Recombinantly expressed in mammalian expression system.
Reactivity	Human
Preservatives	PBS
Form/Appearance	Lyophilized/ Liquid
Storage	Store in -20°C
Tag	C-terminus 6× Histag
Note	For research use only
Application notes	Covid related research
Protein Sequence	V16-Arg-685
Purity	90%, SDS-PAGE (Coomassie blue staining)
Source	HEK
MW	76.5 kDa.
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.