


Product Datasheet

NgR3 Rabbit Monoclonal Antibody (orb1147358)

Description	NgR3 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human NgR3
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	49 kDa
Uniprot ID	Q86UN2
Dilution Range	WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of NgR3 expression...