


Product Datasheet

AdSS 2 Rabbit Monoclonal Antibody (orb1147360)

Description	AdSS 2 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, WB
Immunogen	A synthesized peptide derived from human AdSS 2
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	45 kDa
Uniprot ID	P30520
Dilution Range	WB 1:500-1:2000ICC/IF 1:50-1:200FC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of AdSS 2 expressi...