


Product Datasheet

SPINK5 / LEKTI Rabbit Monoclonal Antibody (orb1147403)

Description	SPINK5 / LEKTI Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human SPINK5 / LEKTI
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	121 kDa
Uniprot ID	Q9NQ38
Dilution Range	WB 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of SPINK5 / LEKTI ...