

Product Datasheet

Thrombomodulin / CD141 Rabbit Monoclonal Antibody (orb1147489)

Description Thrombomodulin / CD141 Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IHC, WB

Immunogen A synthesized peptide derived from human Thrombomodulin / CD141

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Monoclonal

MW 105 kDa

Uniprot ID **P15306**

Dilution Range WB 1:500-1:2000 IHC 1:50-1:200 FC 1:50

Expiration Date 12 months from date of receipt.



Western blot analysis of Thrombomodulin
...