


Product Datasheet

ICAM2 / CD102 Rabbit Monoclonal Antibody (orb1147490)

Description	ICAM2 / CD102 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	A synthesized peptide derived from human ICAM2 / CD102
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	60 kDa
Uniprot ID	P35330
Dilution Range	WB 1:500-1:2000IHC 1:50-1:200IP 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of ICAM2 / CD102 e...