



Product Datasheet

human peroxidasin homolog isoform X1 Antibody (orb1147623)

Description	human peroxidasin homolog isoform X1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Form/Appearance	Liquid
Concentration	500 µg/ml
Storage	At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q92626
Dilution Range	Western blot, 0.5 µg/ml, Hu Immunohistochemistry(Paraffin-embedded Section), 1 µg/ml, Hu
Expiration Date	12 months from date of receipt.



IHC analysis
of
Peroxidasin
homolog
isof...



Western blot
analysis of
Peroxidasin
hom...