



Product Datasheet

Il2ra Antibody (orb1147770)


Description	Il2ra Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC
Immunogen	E.coli-derived mouse Il2ra recombinant protein (Position: E22-I268).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	P01590
Dilution Range	Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Mouse, Rat Immunocytochemistry/Immunofluorescence, 5 µg/ml, Mouse Direct ELISA, 0.1-0.5 µg/ml, Mouse
Expiration Date	12 months from date of receipt.



IHC analysis of Il2ra using anti-Il2ra a...



IHC analysis of Il2ra using anti-Il2ra a...



IF analysis of Il2ra using anti-Il2ra an...