


Product Datasheet

CD118 Antibody (orb1147811)

Description	CD118 (LIFR alpha) is a type I transmembrane glycoprotein of the cytokine receptor family, which ass...
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	human recombinant soluble CD118-Fc
Target	CD118
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P42702
Entrez	3977
Dilution Range	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of JAR cells stained using an...