

Product Datasheet

mCherry Tag Mouse Monoclonal Antibody (9D3) (orb1147836)

Description

Featured Anti-mCherry Tag Mouse

designed for your tagged protein...

Species/Host

Mouse

Conjugation

Unconjugated

Tested

IF, WB

Applications

Immunogen

Recombinant Protein

Preservatives

Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Form/Appearance

Liquid solution

Storage

Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note

For research use only

Application notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:5000), IF (1:200).

Isotype

Mouse IgG

Clonality

Monoclonal

Hazard

Information

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

Expiration Date

12 months from date of receipt.