

Product Datasheet

SuperKine Enhanced Antibody Dilution Buffer (orb1147871)

Description	Antibody Dilution Buffer
Conjugation	Unconjugated
Storage	Stable for one year at 4-8°C from date of shipment. Exposure to the sun or any other intense light can harm the working solution. Short-term exposure to typical laboratory lighting will not harm the working solution.
Note	For research use only
Application notes	SuperKine Enhanced Antibody Dilution Buffer provides an enhanced reagent, which can be directly used to dilute the primary or secondary antibody in an appropriate ratio with this product.
Hazard Information	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.
Expiration Date	12 months from date of receipt.