

Product Datasheet

Thyroxine 4 antibody (orb11479)

Description	Thyroxine 4 antibody
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	BSA conjugated T4
Target	T4/L-Thyroxine
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	affinity purified by Protein A
MW	0.77687kDa
CAS Number	51-48-9
Expiration Date	12 months from date of receipt.