

Product Datasheet

Anti-V5 Antibody (orb1148018)

Description	Anti-V5 Goat Polyclonal Antibody
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Synthetic peptide corresponding to aa 95-108 (GKIPNPPLGLDST) of RNA polymerase alpha subunit of simian virus 5 conjugated to KLH. Host Species: Goat Specificity: Antibody was affinity-purified on the peptide.
Target	V5
Preservatives	0.09% Sodium Azide. PBS, pH 7.2, 0.09% sodium azide.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for at least one (1) year at 2-8C.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Dilution Range	Immunoblotting: use at a dilution of 1:1,000-1:20,000. Immunocytochemistry: use at a dilution of 1:100-1:400. Immunoprecipitation: use 1-4ug/mg lysate. ELISA: for detection use at a dilution of 1:1,000-1:30,000; for coating use at a dilution of 1:100-1:500. These are recommended dilutions. Endusers should determine optimal dilutions for their applications.
Expiration Date	12 months from date of receipt.