

Product Datasheet

Anti-Rabbit IgG H&L Chain Antibody, HRP conjugate (orb1148142)

Description

Anti-Rabbit IgG, H&L Chain,

Antibody

Species/Host

Goat

Reactivity

Rabbit

Conjugation

HRP

Tested

ELISA, ICC, IHC, WB

Applications
Immunogen

Rabbit IgG Host Species: Goat
 Stabilizers: 0.2% Bovine Serum
 Albumin (BSA) Preservatives: 0.1%
 Pro-Clean 400 Buffer: Phosphate-
 buffered saline (PBS) Purification:
 Immunoaffinity chromatography
 Conjugation: Molar enzyme:antibody
 ratio is 4:1

Target

Rabbit IgG H&L Chain, HRP conjugate

Preservatives

0.1% Pro-Clean 400. 0.2% Bovine
 Serum Albumin (rAlbumin), 0.1% Pro-
 Clean 400, Phosphate-buffered saline
 (PBS)

Form/Appearance

Liquid

Concentration

1.0 mg/mL

Storage

Store at 2-8C. Shelf life is one (1)
 year from date of receipt.

Note

For research use only

Clonality

Polyclonal

Purity

Purified

Dilution Range

ELISA: 1:10,000-1:100,000 Western
 Blot: 1:1,000-1:30,000
 Immunohistochemistry: 1:200-1:500
 Immunocytochemistry: 1:200-1:500
 Prepare working dilutions
 immediately before use. These are
 recommended dilutions. User should
 determine optimal dilutions for their
 applications.

Expiration Date

12 months from date of receipt.