

Product Datasheet

Anti-Rabbit IgG H&L Chain Antibody, Alkaline Phosphatase conjugate (orb1148146)

Description

Anti-Rabbit IgG, H&L Chain

Antibody
Species/Host Goat

Reactivity Rabbit

Conjugation AP

Tested Applications ELISA, ICC, IHC, WB

Immunogen Rabbit IgG Host Species: Goat
Stabilizers: 0.2% Bovine Serum Albumin (BSA) Preservatives: 0.09% Sodium Azide Buffer: 50mM HEPES, pH 7.1, 0.1M NaCl, 1mM MgCl₂, 0.1mM ZnCl₂ Purification: Immunoaffinity chromatography Conjugation: Molar enzyme: antibody ratio is 1:1

Target Rabbit IgG H&L Chain, Alkaline Phosphatase conjugate

Preservatives 0.09% Sodium Azide. 0.2% Bovine Serum Albumin (rAlbumin), 0.09% Sodium Azide, 50mM HEPES, pH 7.1, 0.1M NaCl, 1mM MgCl₂, 0.1mM ZnCl₂
Form/Appearance Liquid

Concentration 1.0 mg/mL

Storage Store at 2-8C. Expires 1 year from date of receipt.

Note For research use only

Clonality Polyclonal

Purity Purified

Dilution Range ELISA: 1:1,000-1:10,000 Western Blot: 1:1,000-1:10,000 Immunohistochemistry: 1:200-1:1,000 Immunocytochemistry: 1:200-1:1,000 Prepare working dilutions immediately before use. These are recommended dilutions. User should determine optimal dilutions for their applications.

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558