

Product Datasheet

Anti-Rabbit IgG Fc specific Antibody, Alkaline Phosphatase conjugate (orb1148150)

Description

Anti-Rabbit IgG Fc Fragment, Alk

Species/Host	Goat
Reactivity	Rabbit
Conjugation	AP
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Rabbit IgG Fc fragment Host Species: Goat Stabilizers: 0.2% Bovine Serum Albumin (BSA) Preservatives: 0.09% Sodium Azide Buffer: 50mM HEPES pH 7.1, 0.1M NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂ Purification: Immunoaffinity chromatography Conjugation: Molar enzyme:antibody ratio = 1:1
Target	Rabbit IgG Fc specific, Alkaline Phosphatase conjugate
Preservatives	0.09% Sodium Azide. 0.2% Bovine Serum Albumin (rAlbumin), 0.09% Sodium Azide, 50mM HEPES pH 7.1, 0.1M NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	Store at 2-8C. Shelf life is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Dilution Range	ELISA: 1:5,000-1:50,000 Western Blot: 1:5,000-1:50,000 Immunohistochemistry: 1:100-1:1,000 Immunocytochemistry: 1:200-1:1,000 Prepare working dilutions immediately before use. These are recommended dilutions. User should determine optimal dilutions for their applications.
Expiration Date	12 months from date of receipt.