

Product Datasheet

Anti-Procaïnamide Antibody (orb1148167)

Description	Anti-Procaïnamide Mouse Monoclonal Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Procaïnamide-KLH.
Target	Procaïnamide
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
CAS Number	51-06-9
Dilution Range	These antibodies have been qualified for use in ELISA to detect procaïnamide.
Expiration Date	12 months from date of receipt.