

Product Datasheet

Anti-Mycoplasma bovis Antibody (orb1148227)

Description	Anti-Mycoplasma bovis Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Mycoplasma bovis cells
Target	Mycoplasma bovis
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	These antibodies are stable for at least one (1) year at -20 to -70C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Dilution Range	These antibodies have been qualified for use in ELISA to detect M. bovis. Endusers should determine optimal concentrations for their applications.
Expiration Date	12 months from date of receipt.