

Product Datasheet

Anti-Mycoplasma bovis Mouse Monoclonal Antibody - Library Pack (orb1148228)

Description

Anti-Mycoplasma bovis Mouse

Species/Host

Mouse

Reactivity

Bacteria

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

Mycoplasma bovis cells (ATCC #25025). Host Species: Mouse
Specificity: These antibodies recognize Mycoplasma bovis. They do not cross-react with Mycoplasma bovis genitalium, Mycoplasma californicum, Strep. agalactiae, Staph. aureus, or E. coli.

Target

Mycoplasma bovis

Preservatives

PBS, pH 7.4

Form/Appearance

Liquid

Storage

This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.

Note

For research use only

Clonality

Monoclonal

Purity

Purified

Dilution Range

These antibodies have been qualified for use in ELISA to detect M. bovis. Endusers should determine optimal concentrations for their applications.

Expiration Date

12 months from date of receipt.