

Product Datasheet

Anti-JunB (Phospho-Ser259) Antibody (orb1148322)

Description

Anti-JunB (Phospho-Ser259) Rabbit

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Peptide sequence that includes phosphorylation site of serine 259 (P-B-S(p)-P-I) derived from human JunB and conjugated to KLH. Accession no.: P17275, NP_002220.1 Host Species: Rabbit Specificity: This antibody detects endogenous human, mouse, and rat JunB only when phosphorylated at serine 259.
Target	JunB (Phospho-Ser259)
Preservatives	0.02% Sodium Azide. PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	These antibodies are stable for at least one (1) year at -20C to -70C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P17275
Entrez	165161
Dilution Range	Immunoblotting: use at dilution of 1:500-1:1,000. A band of ~43kDa is detected. Immunohistochemistry: use at dilution of 1:50- 1:100. These are recommended working dilutions. Enduser should determine optimal dilutions for their applications.
Expiration Date	12 months from date of receipt.