

Product Datasheet

Anti-Human IgG H&L Chain Antibody, HRP conjugate (orb1148341)

Description

Anti-Human IgG, H&L Specific-

Species/Host	Goat
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Human IgG Host Species: Goat Stabilizers: 0.2% BSA Preservatives: None Buffer: Phosphate-buffered saline (PBS) containing 0.1% Pro- Clean 400 Purification: Immunoaffinity chromatography Conjugation: Molar enzyme:antibody ratio is 4:1
Target	Human IgG H&L Chain, HRP conjugate
Preservatives	0.1% Pro-Clean 400. 0.2% rAlbumin, None, Phosphate-buffered saline (PBS) containing 0.1% Pro-Clean 400
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Dilution Range	ELISA: 1:10,000-1:100,000 Western Blot: 1:1,000-1:30,000 Immunohistochemistry: 1:200-1:500 Immunocytochemistry: 1:200-1:500 Prepare working dilutions immediately before use. These are recommended dilutions. User should determine optimal dilutions for their applications.
Expiration Date	12 months from date of receipt.