

Product Datasheet

Anti-Human IgA Alpha Chain specific Antibody, Biotin conjugate (orb1148354)

Description

Anti-Human IgA, Alpha Chain

Antibody
Species/Host Goat

Reactivity Human

Conjugation Biotin

Tested Applications ELISA, WB

Immunogen Human IgA Host Species: Goat
Stabilizers: 0.2% Bovine Serum Albumin (BSA) Preservatives: 0.09% Sodium Azide Buffer: Phosphate-buffered saline (PBS) Purification: Immunoaffinity chromatography Conjugation: Biotin

Target Human IgA Alpha Chain specific, Biotin conjugate

Preservatives 0.09% Sodium Azide. 0.2% Bovine Serum Albumin (rAlbumin), 0.09% Sodium Azide, Phosphate-buffered saline (PBS)

Form/Appearance Liquid

Concentration 1.0 mg/mL

Storage Store at 2-8C. Shelf life is one (1) year from date of receipt.

Note For research use only

Clonality Polyclonal

Purity Purified

Dilution Range ELISA: 1:10,000-1:100,000 with enzyme-conjugated streptavidin Western Blot: 1:10,000-1:100,000 with enzyme-conjugated streptavidin Prepare working dilutions immediately before use. These are recommended dilutions. User should determine optimal dilutions for their applications.

Expiration Date 12 months from date of receipt.