

Product Datasheet

Anti-Hemoglobin Antibody (orb1148416)

Description	Anti-Hemoglobin (Human) Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified human hemoglobin
Target	Hemoglobin
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Dilution Range	These antibodies have been qualified for use in ELISA to detect human hemoglobin.
Expiration Date	12 months from date of receipt.