

Product Datasheet

Anti-Haloperidol Antibody (orb1148420)

Description	Anti-Haloperidol Mouse Monoclonal Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Haloperidol conjugated through the tertiary alcohol at the 4 th position of the piperidine ring with a succinic acid linkage to keyhole limpet hemocyanin (KLH). Host Species: Mouse Antibody Class: IgG1 Stabilizers: None Preservatives: None. Available on request.
Target	Haloperidol
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at 2-8C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
CAS Number	52-86-8
Dilution Range	This antibody may be used in immunoassays for detection of haloperidol. Other applications are under investigation.
Expiration Date	12 months from date of receipt.