

Product Datasheet

Anti-HA Antibody, HRP conjugate (orb1148421)

Description	Anti-HA HRP Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Hamster
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	HA (YPYDVPDYA) conjugated to KLH. Host Species: Rabbit Specificity: Antibody was affinity-purified on the peptide.
Target	HA, HRP conjugate
Preservatives	0.1% Pro-Clean 400. PBS, pH 7.2, 0.2% rAlbumin, 0.1% Pro-Clean 400. Molar enzyme/antibody ratio = 4:1.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for at least one (1) year at 2-8C.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Dilution Range	Immunoblotting: use at a dilution of 1:1,000- 1:30,000. Immunocytochemistry: use at a dilution of 1:200-1:500. ELISA: use at a dilution of 1:10,000- 1:100,000. These are recommended dilutions. Endusers should determine optimal dilutions for their applications.
Expiration Date	12 months from date of receipt.