

Product Datasheet

Anti-Glucose Oxidase Mouse Monoclonal Antibody - Library Pack (orb1148437)

Description
Anti-Glucose Oxidase Mouse

Pack - Includes 100 ug of each of the followi...

Species/Host

Mouse

Reactivity

Fungi

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

Purified glucose oxidase from *Aspergillus niger*. Host Species: Mouse Specificity: These antibodies recognize glucose oxidase of *Aspergillus niger*.

Target

Glucose Oxidase

Preservatives

PBS, pH 7.4.

Form/Appearance

Liquid

Storage

This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.

Note

For research use only

Clonality

Monoclonal

Purity

Purified

CAS Number

9001-37-0

Dilution Range

These antibodies have been qualified for use in ELISA to detect glucose oxidase.

Expiration Date

12 months from date of receipt.