

Product Datasheet

Anti-Gentamicin Mouse Monoclonal Antibody - Library Pack (orb1148441)

Description	Anti-Gentamicin Mouse Monoclonal Antibody Library Pack - Includes 100 ug of each of the following an...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Gentamicin conjugated to KLH. Host Species: Mouse Specificity: These antibodies recognize gentamicin.
Target	Gentamicin
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Can be stored at 4C for short-term use.
Note	For research use only
Clonality	Monoclonal
Purity	Purified
CAS Number	1403-66-3
Dilution Range	These antibodies have been qualified for use in ELISA to detect and gentamicin.
Expiration Date	12 months from date of receipt.