

Product Datasheet

Anti-Thyroglobulin Mouse Monoclonal Antibody - Library Pack (orb1148447)

Description
Anti-Thyroglobulin Mouse

- Includes 100 ug of each of the following...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

Purified human thyroglobulin.
Host Species: Mouse Specificity:
12102 is T4-specific and 12104 is
T3-specific. 12101 and 12103
recognize both T3 and T4.

Target

Thyroglobulin

Preservatives

PBS, pH 7.4.

Form/Appearance

Liquid

Storage

This antibody is stable for 1
month at 4°C or one year at -
20°C. Avoid multiple freeze-thaw
cycles.

Note

For research use only

Clonality

Monoclonal

Purity

Purified

Uniprot ID

P01266

Entrez

188450

Dilution Range

These antibodies may be used in
ELISA to detect and quantitate
thyroglobulin.

Expiration Date

12 months from date of receipt.