

Product Datasheet

Anti-Feline Immunoglobulin Antibody (orb1148454)

Description	Anti-Feline Immunoglobulin Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Feline
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Pool of purified feline immunoglobulin from normal sera
Target	Feline Immunoglobulin
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This product is stable for at least one (1) year at -20C to -70C. Reconstituted product should be stored in appropriate aliquots to avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Dilution Range	These antibodies have been qualified for use in ELISA to detect feline immunoglobulin.
Expiration Date	12 months from date of receipt.