

Product Datasheet

Anti-Calnexin Antibody (orb1148462)

Description

Anti-Calnexin Mouse Monoclonal Antibody

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Dictyostelium discoideum AX2 recombinant fragment. Accession no. P27824. Host species: Mouse. Antibody Class: IgG1. Preservative: None.
Target	Calnexin
Preservatives	Lyophilized, 0.1M Tris, 0.1M glycine, 2% sucrose, 1mg/ml.
Form/Appearance	Lyophilized
Concentration	1.0 mg/mL
Storage	This product is stable for at least one (1) year at -20C to -70C. Reconstituted product should be stored in appropriate aliquots to avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P27824
Dilution Range	Immunoblotting: use at 1-10ug/ml. Immunohistochemistry: use at 1-10ug/ml. These are recommended concentrations. Enduser should determine optimal concentrations for their applications.
Expiration Date	12 months from date of receipt.