

Product Datasheet

Anti-CD73/5'Nucleotidase Antibody (orb1148463)

Description

Anti-CD73/5'-Nucleotidase Mouse

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant 5'-nucleotidase (5'-NT). Accession no. P21589. Host species: Mouse. Antibody Class: IgG1. Preservative: None.
Target	CD73/5'-Nucleotidase
Preservatives	Lyophilized, 0.1M Tris, 0.1M glycine, 2% sucrose, 1mg/ml.
Form/Appearance	Lyophilized
Concentration	1.0 mg/mL
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P21589
Dilution Range	ELISA: use at 1-10ug/ml. Immunoblotting: use at 1-5ug/ml. A band of ~63kDa is detected. Immunohistochemistry: use at 1-5ug/ml. These are recommended concentrations. Enduser should determine optimal concentrations for their applications.
Expiration Date	12 months from date of receipt.