

## Product Datasheet

### Anti-Feline Immunoglobulin Antibody (orb1148464)

|                            |                                                                                                                                                    |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Anti-Feline Immunoglobulin Mouse<br>Monoclonal Antibody                                                                                            |
| <b>Species/Host</b>        | Mouse                                                                                                                                              |
| <b>Reactivity</b>          | Feline                                                                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                       |
| <b>Tested Applications</b> | ELISA                                                                                                                                              |
| <b>Immunogen</b>           | Pool of purified feline immunoglobulin from normal sera                                                                                            |
| <b>Target</b>              | Feline Immunoglobulin                                                                                                                              |
| <b>Preservatives</b>       | PBS, pH 7.4.                                                                                                                                       |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                             |
| <b>Storage</b>             | These antibodies are stable for at least one (1) year at -20C to -70C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                              |
| <b>Isotype</b>             | IgG1                                                                                                                                               |
| <b>Clonality</b>           | Monoclonal                                                                                                                                         |
| <b>Purity</b>              | Purified                                                                                                                                           |
| <b>Dilution Range</b>      | These antibodies have been qualified for use in ELISA to detect feline immunoglobulin.                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                    |