

Product Datasheet

Anti-FAK (Phospho-Tyr861) Antibody (orb1148470)

Description

Anti-FAK (Phospho-Tyr861) Rabbit

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide sequence that includes phosphorylation site of tyrosine 861 (H-I-Y(p)-Q-P) derived from human FAK. Host Species: Rabbit Accession no.: Q05397, NP_005598.3 Specificity: This antibody detects endogenous human FAK only when phosphorylated at tyrosine 861.
Target	FAK (Phospho-Tyr861)
Preservatives	0.02% Sodium Azide. PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	These antibodies are stable for at least one (1) year at -20C to -70C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Dilution Range	Immunoblotting: use at dilution of 1:500-1:1,000. A band of ~125kDa is detected. These are recommended working dilutions. Enduser should determine optimal dilutions for their applications.
Expiration Date	12 months from date of receipt.