

Product Datasheet

Anti-Tau Antibody (orb1148537)

Description	Anti-Tau Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide derived from the full-length sequence of human Tau
Target	Tau
Preservatives	0.09% Sodium Azide. PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P10636
Dilution Range	ELISA: use at 1-10ug/ml with human Tau on the solid phase. Immunoblotting: use at 1-5ug/ml. A band of ~45.8kDa is detected.
Expiration Date	12 months from date of receipt.