

Product Datasheet

Anti-Ankyrin-B Antibody (orb1148590)

Description

Anti-Ankyrin-B Mouse Monoclonal Antibody

Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide corresponding to aa 203-496 of human Ankyrin-B (accession no. NP_001120965.1) Host Species: Mouse Antibody Class: IgG1
Target	Ankyrin-B
Preservatives	0.09% Sodium Azide. PBS, pH 7.4, 50% glycerol, 0.09% sodium azide. Purified by Protein G affinity chromatography.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This product is stable for at least 1 year at -20C. Freeze in multiple aliquots to avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q01484
Dilution Range	Immunoblotting: use at 1-2ug/ml. A band of >200kDa is detected. Immunohistochemistry: use at 1-5ug/ml. These are recommended concentrations. Enduser should determine optimal concentrations for their applications. Positive control: rat brain membranes.
Expiration Date	12 months from date of receipt.