

Product Datasheet

Anti-Bassoon Antibody (orb1148598)

Description	Anti-Bassoon Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Fusion protein corresponding to aa786-1041 (N-terminal) of human Bassoon (accession no. NM_003458.3). Host Species: Rabbit Antibody Class: Polyclonal
Target	Bassoon
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q9UPA5
Dilution Range	Immunoblotting: use at 1ug/ml. Immunohistochemistry: use at 1-10ug/ml. Immunofluorescence: use at 1-10ug/ml. Positive control: Rat brain tissue lysate These are recommended concentrations; enduser should determine optimal concentrations for their applications.
Expiration Date	12 months from date of receipt.