

Product Datasheet

Anti-PSD95 Antibody (orb1148613)

Description	Anti-PSD95 Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Bovine, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant rat PSD95 Host Species: Mouse Gene ID: 29495 Accession No.: NP_062567.1, P31016 Specificity: This antibody recognizes mouse, rat, pig, and human PSD95.
Target	PSD95
Preservatives	0.09% Sodium Azide. PBS, pH 7.4, 50% glycerol, 0.09% sodium azide.
Form/Appearance	Liquid
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P31016
Dilution Range	Immunoblotting: use at 1ug/ml. A band of ~100 kDa is detected. Immunocytochemistry: use at 1-10ug/ml. Immunofluorescence: use at 10ug/ml. These are recommended concentrations. Endusers should determine optimal concentrations for their application.
Expiration Date	12 months from date of receipt.