

Product Datasheet

Anti-Rab5 Antibody (orb1148620)

Description	Anti-Rab5 Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Host Species: Rabbit Human Rab5 synthetic peptide (identical to canine Rab5) conjugated to KLH Antibody Class: Polyclonal Preservatives: None.
Target	Rab5
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid repeated freezing and thawing.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q6FI44
Dilution Range	Immunoblotting: use at 1ug/ml. A band of ~24-26 kDa is detected. Positive control: HeLa cell lysate.
Expiration Date	12 months from date of receipt.