

Product Datasheet

Anti-Calnexin (NT) Antibody (orb1148621)

Description	Anti-Calnexin Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Host Species: Rabbit A 19 amino acid synthetic peptide corresponding to canine calnexin conjugated to KLH. Antibody Class: Polyclonal Preservatives: None.
Target	Calnexin (NT)
Preservatives	Whole Antiserum
Form/Appearance	Liquid
Storage	Stable for one year at -20°C. Store in appropriate aliquots to avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Serum
Uniprot ID	P27824
Dilution Range	Immunoblotting: use at 1:5,000-1:10,000 dilution. A band of ~90 kDa is detected. These are recommended concentrations. User should determine optimal concentrations for their application. Positive control: HeLa cell lysate.
Expiration Date	12 months from date of receipt.