

Product Datasheet

Anti-IRAK-4 Antibody (orb1148622)

Description

Anti-IRAK-4 Rabbit Polyclonal Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

A synthetic peptide corresponding to amino acids at the carboxy terminus of human IRAK-4. Host Species: Rabbit Preservative: 0.02% sodium azide

Target

IRAK-4

Preservatives

PBS, pH 7.4.

Form/Appearance

Liquid

Concentration

1.0 mg/mL

Storage

This product does not contain preservative. Lyophilized antibody is stable at 4 C for 2 years. Store the reconstituted solution (2.0mg/ml) in aliquots at -20 C for one year or at 4 C for 6 months after adding 0.1% sodium azide. Dilutions of the stock solution should not be stored. Avoid repeated freeze-thaw cycles.

Note

For research use only

Clonality

Polyclonal

Purity

Purified

Uniprot ID
Q9NWZ3
Dilution Range

Immunoblotting: Use at 1-4ug/ml. A band of approximately 50kDa is detected.
Immunocytochemistry: Use at 10ug/ml.
Immunofluorescence: Use at 10ug/ml. These are recommended concentrations. Enduser should determine optimal concentrations for their applications.

Expiration Date

12 months from date of receipt.