

Product Datasheet

Anti-Calpastatin Antibody (orb1148623)

Description	Anti-Human Calpastatin Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant human muscle- type calpastatin. Host Species: Mouse Clone Number: CSL1-5 Subtype: IgG1 Stabilizers: 1% BSA Preservatives: None
Target	Calpastatin
Preservatives	10 mM PBS, pH 7.4, 1% rAlbumin, lyophilized
Form/Appearance	Lyophilized
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P20810
Dilution Range	ELISA: 1-10ug/ml with immobilized antigen. Western Blot: 1-10ug/ml, denaturing and non-denaturing conditions. Immunohistochemistry: 1-10ug/ml, frozen tissue sections.
Expiration Date	12 months from date of receipt.