

Product Datasheet

Anti-DENTT Antibody (orb1148624)

Description	Anti-DENTT Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human DENTT protein. Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide
Target	DENTT
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q9H2G4
Dilution Range	Immunoblotting: use at 1:500-1:1,000 dilution. In immunoblots, a band of approximately 77 kD is detected.
Expiration Date	12 months from date of receipt.