

Product Datasheet

Anti-GABABR2 Antibody (orb1148627)

Description	Anti-GABA(B)R2 Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide corresponding to aa 42-54 (TRGAPRPPSSPP) of the rat GABABR2 subunit. Host Species: Rabbit
Target	GABABR2
Preservatives	0.02% Sodium Azide. PBS, pH 7.4, 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	O75899
Dilution Range	Immunoblotting: use at 1:250 – 1:1,000 dilution. A band of ~105kDa is detected.
Expiration Date	12 months from date of receipt.