

Product Datasheet

Anti-GABABR1 Antibody (orb1148628)

Description	Anti-GABA(B)R1 Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to aa 33-54 of the GABABR1 subunit. Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide
Target	GABABR1
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q9UBS5
Dilution Range	Immunoblotting: use at 1:500-1:1,000 dilution. In immunoblots, a band of approximately 120 kD is detected.
Expiration Date	12 months from date of receipt.