

Product Datasheet

Anti-GABA Receptor Associated Protein Antibody (orb1148629)

Description	Anti-GABARAP Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide derived from the N-terminus of GABARAP. Host Species: Rabbit Preservatives: 0.02% sodium azide.
Target	GABA Receptor Associated Protein
Preservatives	PBS.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid repeated freezing and thawing.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	O95166
Dilution Range	Immunoblotting: use at 1:500-1:1,000 dilution. ELISA: use at 1:1,000-1:5,000 dilution.
Expiration Date	12 months from date of receipt.